

Sure Puncture Needle

(Initial Puncture Needle - Three Part)



Type
Three Part

Size (in G)
21/16

Length (in cm)
20

Intended Use:

Intended to make the initial percutaneous puncture to access the kidney and to enable the introduction of a guidewire during urological procedures (PCN, PCNL, Mini PERC & Micro PERC).

Key Features:

- Siliconised needle for less trauma
- Echogenic tip for better visualization during ultrasound
- 1 cm marking

Advantages of Three Part IP Needle:

- Combination of 21 G Bevel tip and 16 G Cannula Needle
- Inner 21 G provides minimal site puncture and contrast injection to ensure desired puncture location
- Coaxial needle provides secured access
- Can introduce working guidewire with 21 G puncture