

## Sure Puncture Needle (Initial Puncture Needle Three Part)

### Intended Use:

Intended to make the initial percutaneous puncture to create access and to enable the introduction of guidewire during percutaneous interventional procedures



Type	Size (in G)	Length (in cm)
Three part	21/16	20

### Key Features and Benefits:

- **Siliconized needle:** Less trauma during insertion
- **Echogenic tip:** Better visualization during ultrasound guided procedures
- **1 cm marking:** Helps in gauging depth during insertion
- **Combination of 21 G Bevel tip and 16 G Cannula Needle:** Provides optimal sharpness for initial puncture along with a larger cannula for guidewire introduction
- **Compatible with 21 G puncture:** Enables the introduction of a working guidewire
- **Inner 21 G:** Minimal site puncture and allows contrast injection to ensure desired puncture location
- **Coaxial needle:** Provides secured access during the procedure