

# SURE PUNCTURE<sup>™</sup> NEEDLE

## (INITIAL PUNCTURE NEEDLE - THREE PART)

Intended to make the initial percutaneous puncture to create access and to enable the introduction of a guidewire during percutaneous interventional procedures.

### Key Features

- Siliconised needle for less trauma
- Echogenic tip for better visualization during ultrasound
- 1 cm marking



### Type

Three Part

### Size (in G)

21/16

### Length (in cm)

20

### Advantage of Sure Puncture<sup>™</sup> Needle

- Combination of 21 G Bevel tip and 16 G Cannula Needle
- Inner 21 G provides minimal site puncture & contrast injection to ensure desired puncture location
- Coaxial needle provides secured access
- Can introduce working guidewire with 21 G puncture