


Accura Kit


Percutaneous Access Set

Used for single puncture percutaneous access to facilitate placement of 0.038" Guide wire for interventional radiology procedures.



Key Features & Benefits

- **Non Kinkable Introducer Sheath:** Prevents bending or collapse during insertion for smoother navigation and higher procedural success
- **High Radiopaque Polymer Band:** Ensures excellent visibility under fluoroscopy or X-ray, enabling precise placement even in complex anatomy
- **Cope Mandril Guidewire with Floppy Tip:** Minimizes trauma and enhances safety in delicate structures
- **21G Needle with 15 cm Length:** Fine access needle enables minimal trauma during initial puncture
- **Hydrophilic Coating:** For smooth insertion



Set Includes

- 21G Chiba needle with Echogenic tip
- 0.018" X 60 cm Cope Mandril Wire
- Introducer Sheath with Dilator

Specifications

Size (in Fr)
4, 5 & 6

Length (in cm)
19


360° Product Experience



Accura Micro Kit




Percutaneous Access Set

Used for single puncture percutaneous access to facilitate placement of 0.038" Guide wire for interventional radiology procedures.



Key Features & Benefits

- **Non Kinkable Introducer Sheath:** Prevents bending or collapse during insertion for smoother navigation and higher procedural success
- **High Radiopacity:** Ensures excellent visibility under fluoroscopy or X-ray, enabling precise placement even in complex anatomy
- **Cope Mandril Guidewire with Floppy Tip:** Minimizes trauma and enhances safety in delicate structures
- **21G Needle with 7 cm Length:** Fine access needle enables minimal trauma during initial puncture



Set Includes

- 21G Introducer Needle
- 0.018" X 40 cm Cope Mandril Wire
- Introducer Hub
- Introducer Sheath with Dilator

Specifications

Size (in Fr)
4 & 5

Length (in cm)
10

BLUENEEM[®]
Clinicians Trust. Patients Wellbeing

Blue Neem Medical Devices Pvt Ltd
Plot Nos 270 & 271, Road No 5
Harohalli Industrial Area, Phase II
Kanakapura Taluk, Ramanagara - 562112
Karnataka. INDIA.

📞 +91 97399 72855

✉️ enquiry@blueneem.com

🌐 www.blueneem.com