

# ACCCURA KIT

## PERCUTANEOUS ACCESS SET

Used for single puncture percutaneous access to facilitate placement of 0.038 inch diameter Guidewire for interventional radiology procedures.

### Key Features

- A radiopaque band located at the tip of the sheath helps to identify precise location of the sheath's distal tip for accurate positioning
- The hydrophilic coating on sheath when activated, becomes highly lubricious enabling smooth insertion by reducing friction
- Sheath lumen allows both 0.018 inch cope mandril wire and 0.038 inch wire together
- Set includes a chiba needle with echogenic tip which enhances visibility under ultrasound



### Set Includes

6 Fr introducer sheath, 21 G x 15 cm chiba needle with echogenic tip and 0.018 inch nitinol cope mandril wire