

## Amplatz Renal Dilator (Renal Dilator With Amplatz Sheath)



**Dilator Size**  
(OD of Dilator in Fr)  
12, 14, 16, 18, 20, 22, 24  
26, 28 & 30

**Length (in cm)**  
30

**Amplatz Sheath Size**  
(ID of Sheath in Fr)  
12, 14, 16, 18, 20, 22, 24  
26, 28 & 30

**Length (in cm)**  
17

### Intended Use:

Used for dilation of a percutaneous tract and setup of a working channel into the kidney for PCNL procedures by placing a suitable amplatz sheath.

### Key Features:

- Gradual taper designed for atraumatic placement and maneuverability
- Smooth surface for ease of introduction
- Optimal radiopacity for better visualization