

AMPLATZ RENAL DILATION SET

Intended for stretching or enlarging the fascial tissue prior to a nephrostomy procedure and to introduce the amplatz sheath which is intended for percutaneous dilation of tissue surrounding a dilator previously placed and used as a working channel during urological procedures.

Set Includes

- 8 Fr guiding catheter
- 10 Fr radiopaque dilation sheath
- Radiopaque fascial dilators of sizes from 6 Fr – 10 Fr of length 20 cm
- Radiopaque renal dilators of sizes from 12 Fr – 30 Fr of length 30 cm
- Teflon amplatz sheath for renal dilators of sizes 24, 26, 28 & 30 Fr of length 17 cm
- Packed in tray

