

AMPLATZ RENAL DILATOR

(RENAL DILATOR WITH AMPLATZ SHEATH)

Used for dilation of a percutaneous tract and setup of a working channel into the kidney for PCNL procedures by placing a suitable amplatz sheath.

Key Features

- Gradual taper designed for atraumatic placement and maneuverability
- Smooth surface for ease of introduction
- Optimal radiopacity for better visualization

Dilator Size (OD of Dilator in Fr)

12, 14, 16, 18, 20, 22, 24, 26, 28 & 30

Length (in cm)

30

Amplatz Sheath Size (ID of Sheath in Fr)

12, 14, 16, 18, 20, 22, 24, 26, 28 & 30

Length (in cm)

17

