


# Denssicath<sup>®</sup>

## Multipurpose Catheter

Used for percutaneous drainage of thick and dense abscesses in a variety of drainage procedures.

### Procedures:

- Liver Abscess
- Biliary
- Pancreatic Pseudocyst
- Cholecystostomy(Gallbladder)
- Abdominal Cavity
- Peritoneal Cavity
- Ascites
- Pleural Effusion
- Pneumothorax
- Haemothorax
- Pericardial Drainage




# Key Features And Benefits



The length markings on catheters help determine the insertion depth



Larger eye holes for efficient drainage of thick abscesses, minimizing the risk of clogging




Hydrophilic coating ensures easy and smooth catheter insertion




Ideal tip stiffness to prevent buckling during insertion




Pigtail tip configuration ensures secure retention of the distal tip



Locking mechanism ensures secure placement and prevents accidental catheter slippage



Malecot Flower design is optimally suited for effective drainage of thick, viscous fluids



Straight catheter facilitates effective external biliary drainage in cases of bile duct obstructions

# Product Range

## Denssicath® Multi-Purpose Drainage Catheter



### Package Includes

- Drainage catheter with hydrophilic coating
- Metal and plastic straightener
- Drainage bag connector


## Denssicath® Multi-Purpose Drainage Catheter with Needle



### Package Includes

- Drainage catheter with hydrophilic coating
- Metal and plastic straightener
- Metal trocar
- Drainage bag connector

## Denssicath® Multi-Purpose Drainage Catheter Set



### Set Includes

- Scalpel blade #11
- Chiba needle 18G
- Amplatz Ultra Stiff® guidewire
- Fascial dilators
- Drainage catheter with hydrophilic coating
- Metal and plastic straightener
- Catheter fixing disc
- Drainage bag Connector

## Product Specifications

Product Name	Size (in Fr)	Length (in cm)	Types
Denssicath® Multi-Purpose Drainage Catheter	6, 7, 8, 9, 10, 12, 14 & 16	22, 25, 30 & 35	Pigtail
	7, 8, 9, 10, 12, 14 & 16	22, 25, 30 & 35	Pigtail with Locking Mechanism
	7, 8, 9, 10, 12, 14 & 16	22, 25, 30 & 35	Pigtail with Needle
	8, 10, 12, 14, 16, 18 & 20	22, 25, 30 & 35	Malecot
	8, 10, 12, 14, 16, 18 & 20	22, 25, 30 & 35	Malecot with Needle
	6, 7, 8, 9, 10, 12, 14 & 16	22, 25, 30 & 35	Straight
	6, 7, 8, 9, 10, 12, 14 & 16	22, 25, 30 & 35	Straight with Needle
Denssicath® Multi-Purpose Drainage Catheter Set	6, 7, 8, 9, 10, 12, 14 & 16	22, 25, 30 & 35	Pigtail
	7, 8, 9, 10, 12, 14 & 16	22, 25, 30 & 35	Pigtail with Locking Mechanism
	8, 10, 12, 14, 16, 18 & 20	22, 25, 30 & 35	Malecot

**BLUENEEM**®

CLINICIANS' TRUSTED PARTNER



**Blue Neem Medical Devices Private Limited**

Plot No 270 & 271, Road No.5  
Harohalli Industrial Area - Phase II  
Kanakapura Taluk, Ramanagara District  
Karnataka(INDIA) - 562 112



+91 97399 72855



enquiry@blueneem.com

[www.blueneem.com](http://www.blueneem.com)

