

Biliary Drainage Catheters

Denssiflo[®]

PTBD Catheter

Used for percutaneous placement, for trans hepatic biliary drainage, up to the duodenum for **internal-external** drainage of bile.




Denssiflo[®]


PTBD Catheter




Available
in three hole
configurations
16, 24 & 32



Optimal Retention



Position Marker for Better Visibility & Placement



Larger Eye Holes for Better Drainage



Hydrophilic Coating for Easy Insertion

For Internal-External Biliary Drainage

Denssiflo® PTBD Catheter Set Standard

Set Contains

- #11 Scalpel Blade
- Chiba Needle 18 G
- Amplatz Ultra Stiff® Guidewire
- Fascial Dilators
- Plastic Straightener
- PTBD Catheter with Position Marker and Hydrophilic Coating
- Catheter Fixing Disc
- Drainage Bag Connector

Types

Standard & Ring Biliary

Size (in Fr)

6, 8, 8.5, 10, 12
14 & 16

Length (in cm)

30, 35, 40
45 & 50



Denssiflo® PTBD Catheter Set Micro Puncture

Set Contains

- #11 Scalpel Blade
- Chiba Needle 21 G
- 0.018 inch Cope Mandril Wire
- 6Fr Coaxial Introducer Sheath
- Amplatz Ultra Stiff® Guidewire
- Fascial Dilators
- Plastic Straightener
- PTBD Catheter with Position Marker and Hydrophilic Coating
- Catheter fixing disc
- Drainage bag connector

Size (in Fr)

6, 8, 8.5, 10, 12
14 & 16

Length (in cm)

30, 35, 40
45 & 50



For External Biliary Drainage


Denssicath® Multipurpose Drainage Catheter - Straight

Size (in Fr)

6, 7, 8, 9, 10, 12
14 & 16

Length (in cm)

22, 25, 30
& 35



Product Specifications

Product Name	Type	Size (in Fr)	Length (in cm)	Hole Configuration
Denssiflo ® PTBD Catheter with Position Marker	Pigtail	6, 8, 8.5, 10, 12, 14 & 16	30, 35, 40, 45 & 50	16 or 32 holes
Denssiflo ® PTBD Catheter with Position Marker	Ring Biliary	8, 8.5, 10, 12, 14 & 16	35, 40, 45 & 50	24 holes
Denssiflo ® PTBD Catheter with Position Marker & Locking Mechanism	Pigtail	8, 10, 12, 14 & 16	35, 40 & 45	16 holes
Denssiflo ® PTBD Catheter Set with Position Marker Standard	Pigtail	6, 8, 8.5, 10, 12, 14 & 16	30, 35, 40, 45 & 50	16 or 32 holes
Denssiflo ® PTBD Catheter Set with Position Marker Micro Puncture	Pigtail	6, 8, 8.5, 10, 12, 14 & 16	30, 35, 40, 45 & 50	16 or 32 holes
Denssicath ® Multipurpose Drainage Catheter	Straight	6, 7, 8, 9, 10, 12, 14 & 16	22, 25, 30 & 35	NA



CLINICIANS' TRUSTED PARTNER



Blue Neem Medical Devices Private Limited

Plot No 270 & 271, Road No.5
Harohalli Industrial Area - Phase II
Kanakapura Taluk, Ramanagara District
Karnataka(INDIA) - 562 112



+91 97399 72855



enquiry@blueneem.com

www.blueneem.com

360° Product Experience

