



Denssima[®]

Centesis Catheter


Used to safely and effectively drain clear ascites and pleural effusions while ensuring secure and precise placement




Key Features And Benefits




Hydrophilic coating for easy and smooth insertion




Soft and safe atraumatic tip, eliminating tissue perforation risk



Multiple side holes ensure faster drainage




Suitable for indwelling upto 30 days, avoiding repeated punctures in frequent drainage cases



The length markings on catheters help determine the insertion depth

Product Range

Denssima[®] Centesis Drainage Set



Set Includes

- Drainage tube with hydrophilic coating
- 18G x 7cm introducer needle
- Introducer guidewire (0.038" x 50cm)
- Drainage bag connector

Size (in Fr)
5 & 6

Length (in cm)
10 & 16

Type
Straight & Pigtail

Denssima[®] Single Step Centesis Drainage Set



Set Includes

- Drainage tube with hydrophilic coating
- Puncture needle with stopper
- Drainage bag connector

Size (in Fr)
5 & 6

Length (in cm)
10 & 16

Type
Straight & Pigtail

Product Specifications

Product Name	Type	Size (in Fr)	Length (in cm)
Denssima® Centesis Drainage Catheter Set	Pigtail	5	16
	Straight	5	16
	Pigtail	6	16
	Straight	6	16
Denssima® Single Step Centesis Drainage Catheter Set	Straight	5	10
	Straight	5	16
	Pigtail	5	16
	Straight	6	10
	Straight	6	16
	Pigtail	6	16

BLUENEEM®
CLINICIANS' TRUSTED PARTNER



Blue Neem Medical Devices Private Limited
Plot No 270 & 271, Road No.5
Harohalli Industrial Area - Phase II
Kanakapura Taluk, Ramanagara District
Karnataka(INDIA) - 562 112



+91 97399 72855



enquiry@blueneem.com

www.blueneem.com

