


PrimacoreTM

Biopsy Gun - Semi Automatic

Designed for **smooth** and **efficient** collection of **exact, high quality** core samples from soft tissues and organs.



Smooth Stylet Insertion

Ensures minimal or no tissue displacement



Echogenic Needle Tip

Excellent tip visibility under ultrasound guidance

Depth Markings

To facilitate the precise placement of the needle at the desired depth

Penetration Depth Indicator

Confirms primed penetration depth of 10mm or 20mm

Ergonomic Handle Design

Ergonomic handle with multiple grip options

Fire-ready Indicator

Designed to reduce the risk of premature firing

SMOOTH access:

- Ergonomic handle with multiple grip options
- Excellent tip geometry to reduce penetration force and tissue displacement
- Exemplary tip visibility under ultrasound guidance

EFFICIENT operation:

- Convenient adjustable throw lengths – 10mm and 20mm
- Smooth stylet insertion
- Siliconized needles to reduce needle deflection
- High Stylet stiffness avoids needle bending

EXACT, High-quality sample collection:

- Penetration depth indicator – confirms primed penetration depth
- Semi-automatic firing – aids visual verification of sample notch placement
- Superlative tactile feedback of tissue penetration and collection

Set Includes:

- Primacore™ Biopsy Instrument
- Compatible Coaxial
- Trocar and Blunt Tip Stylet
- Depth Stopper
- 10 mm Adapter

Product Specifications

Size (G)	Length (cm)
14	10, 16 & 20
16	10, 16 & 20
18	10, 16, 20 & 25
20	10, 16 & 20

