


## Bluthane<sup>®</sup> Ureteral Stent



### Types

- Both Ends Open
- One End Closed
- Both Ends Multiloop
- One End Multiloop
- With or Without Hydrophilic coating
- With or Without String

### Package Includes

Stent with pusher

### Also available in set, which includes

Stent with pusher & guidewire

### Intended Use:

Intended to be implanted for drainage of urine from renal pelvis to urinary bladder.

### Key Features:

- Thermosensitive, advanced PU material for enhanced patient comfort
  - Stent will be stiffer at room temperature promoting pushability for stent deployment
  - Material becomes softer at body temperature
- Optimal loop shape and high memory retention reduces the risk of migration
- Optimal radiopacity for better visualization

### Standard

#### Size (in Fr)

4.5, 4.7, 4.8, 5, 5.5, 6, 6.5, 7 & 8

#### Length (in cm)

12, 14, 16, 18, 20, 22, 24, 26, 28 & 30

### Multiloop

#### Size (in Fr)

4.7, 4.8, 5, 6, 7 & 8

#### Length (in cm)

18-24 & 24-30

### Pediatric

#### Size (in Fr)

2.8, 3, 3.5, 3.7, 3.8 & 4

#### Length (in cm)

6, 8, 10, 12, 14 & 16