

CYTOCORE

BIOPSY GUN - SEMI AUTOMATIC

Used for obtaining histological core samples from soft tissue and organs.

Key Features

- Echogenic and atraumatic needle tip
- Ergonomic and light weight
- Cannula with depth marking
- Easy and one hand operated with two stage arming
- Package includes Coaxial Needle

BIOPSY GUN - SEMI AUTOMATIC

Size	Length
14 G	10 & 16 cm
16 G	10, 16 & 20 cm
18 G	10, 16, 20 & 25 cm
20 G	10, 16 & 20 cm

