

Delta Basket

Intended Use:

Steerable Basket for rigid nephroscope for stone extraction from calices which are inaccessible due to acute angle



Size (in Fr)
7

Length (in cm)
40

Key Features and Benefits:

- **Steerable Design:** Access to inner calyx, including acute angles.
- **Angled or Straight Use:** Versatile stone manipulation.
- **Unique Basket Shape:** Minimizes stone migration.
- **Atraumatic Design:** Prevents infundibular injury and reduces bleeding.
- **Single Tract Access:** Avoids the need for multiple percutaneous nephrostomy tracts
- **Rigid Scope Compatibility:** Eliminates the need for a flexible nephroscope in some cases
- **Clot Removal:** Helps to remove large blood clots