

## Delta Basket



Size (in Fr)

7

Length (in cm)

40

### Intended Use:

Steerable Basket for rigid nephroscope for stone extraction from calices which are inaccessible due to acute angle.

### Design Features:

- Can be angled inside kidney to reach inner calyx
- Can be used for straight path stone manipulation and also can be angled both directions
- Unique basket shape which prevents stone migration

### Procedural Features:

- Prevents injury to the infundibulum and reduce bleeding
- Avoids the need for multiple tracts
- Can be used without the assistance of flexible nephroscope
- Helps to remove big blood clots