


Denssicath[®]

Multipurpose Catheter

Used for percutaneous drainage of thick and dense abscesses in a variety of drainage procedures.

Procedures:

- Liver Abscess
- Biliary
- Pancreatic Pseudocyst
- Cholecystostomy (Gallbladder)
- Abdominal Cavity
- Peritoneal Cavity
- Ascites
- Pleural Effusion
- Pneumothorax
- Haemothorax
- Pericardial Drainage




Key Features And Benefits




The length markings on catheters help determine the insertion depth



Larger eye holes for efficient drainage of thick abscesses, minimizing the risk of clogging




Hydrophilic coating ensures easy and smooth catheter insertion



Ideal tip stiffness to prevent buckling during insertion




Pigtail tip configuration ensures secure retention of the distal tip



Locking mechanism ensures secure placement and prevents accidental catheter slippage



Malecot Flower design is optimally suited for effective drainage of thick, viscous fluids



Straight catheter facilitates effective external biliary drainage in cases of bile duct obstructions

Product Range

Denssicath® Multi-Purpose Drainage Catheter



Package Includes

- Drainage catheter with hydrophilic coating
- Metal and plastic straightener
- Drainage bag connector


Denssicath® Multi-Purpose Drainage Catheter with Needle



Package Includes

- Drainage catheter with hydrophilic coating
- Metal and plastic straightener
- Metal trocar
- Drainage bag connector

Denssicath® Multi-Purpose Drainage Catheter Set



Set Includes

- Scalpel blade #11
- Chiba needle 18G
- Amplatz Ultra Stiff® guidewire
- Fascial dilators
- Drainage catheter with hydrophilic coating
- Metal and plastic straightener
- Catheter fixing disc
- Drainage bag Connector

Product Specifications

Product Name	Size (in Fr)	Length (in cm)	Types
Denssiccath® Multi-Purpose Drainage Catheter	6, 7, 8, 9, 10, 12, 14 & 16	22, 25, 30 & 35	Pigtail
	7, 8, 9, 10, 12, 14 & 16	22, 25, 30 & 35	Pigtail with Locking Mechanism
	7, 8, 9, 10, 12, 14 & 16	22, 25, 30 & 35	Pigtail with Needle
	8, 10, 12, 14, 16, 18 & 20	22, 25, 30 & 35	Malecot
	8, 10, 12, 14, 16, 18 & 20	22, 25, 30 & 35	Malecot with Needle
	6, 7, 8, 9, 10, 12, 14 & 16	22, 25, 30 & 35	Straight
	6, 7, 8, 9, 10, 12, 14 & 16	22, 25, 30 & 35	Straight with Needle
Denssiccath® Multi-Purpose Drainage Catheter Set	6, 7, 8, 9, 10, 12, 14 & 16	22, 25, 30 & 35	Pigtail
	7, 8, 9, 10, 12, 14 & 16	22, 25, 30 & 35	Pigtail with Locking Mechanism
	8, 10, 12, 14, 16, 18 & 20	22, 25, 30 & 35	Malecot



Clinicians Trust. Patients Wellbeing



Blue Neem Medical Devices Private Limited

Plot No 270 & 271, Road No.5
Harohalli Industrial Area - Phase II
Kanakapura Taluk, Ramanagara District
Karnataka(INDIA) - 562 112



+91 97399 72855



enquiry@blueneem.com

www.blueneem.com



360° Product
Experience