


DENSSICATH[®]

MULTI-PURPOSE DRAINAGE CATHETER - MALECOT WITH NEEDLE

Used for percutaneous drainage in a variety of drainage applications (e.g., nephrostomy, biliary and abscess) by direct stick technique.

Key Features

- Integral flower design matches with the catheter size when closed and cause less trauma during insertion
- Trocar design and catheter design prevents buckling of catheter during the procedure
- The Malecot flower design is ideally suitable for the drainage of thick viscous fluids
- Hydrophilic coating for easy insertion



Size (in Fr)

8, 10, 12, 14, 16, 18, 20 & 22

Length (in cm)

20, 25, 30 & 35

Package Includes

Malecot catheter with hydrophilic coating, metal & plastic straightner & drainage bag connector