

DENSSICATH[®]

MULTI-PURPOSE DRAINAGE CATHETER - PIGTAIL

Used for percutaneous drainage in a variety of drainage applications (e.g., nephrostomy, biliary and abscess) by Seldinger access technique.

Key Features

- The pigtail tip configuration helps in the retention of the distal tip in renal pelvis
- Eye holes for better drainage
- Hydrophilic coating for easy insertion



Size (in Fr)

6, 7, 8, 9, 10, 12, 14 & 16

Length (in cm)

22, 25, 30 & 35

Package Includes

Pigtail drainage catheter with hydrophilic coating, metal & plastic straightener and drainage bag connector