

DENSSICATH[®]

MULTI-PURPOSE DRAINAGE CATHETER - PIGTAIL WITH LOCKING MECHANISM

Used for percutaneous drainage in a variety of drainage applications (e.g., nephrostomy, biliary and abscess) by Seldinger access technique. The Locking mechanism prevents accidental slippage of the catheter.

Key Features

- The pigtail tip configuration helps in the retention of the distal tip
- The locking mechanism prevents accidental slippage of the catheter
- Eye holes for better drainage
- Hydrophilic coating for easy insertion



Size (in Fr)

7, 8, 9, 10, 12, 14 & 16

Length (in cm)

22, 25, 30 & 35

Package Includes

Pigtail drainage catheter with locking mechanism & hydrophilic coating, metal & plastic straightener and drainage bag connector.

Note: Once the catheter is deployed, pull the string, wind and lock to the hub of the catheter