

DENSSICATH[®]

MULTI-PURPOSE DRAINAGE CATHETER - PIGTAIL WITH NEEDLE

Used for percutaneous drainage in a variety of drainage applications (e.g., nephrostomy, biliary and abscess) by direct stick technique.

Key Features

- The pigtail tip configuration helps in the retention of the distal tip
- Eye holes for better drainage
- Hydrophilic coating for easy insertion
- Trocar design and catheter design prevents buckling of catheter during the procedure
- Available with or without locking mechanism



Size (in Fr)

7, 8, 9, 10, 12, 14 & 16

Length (in cm)

22, 25, 30 & 35

Package Includes

Pigtail drainage catheter with hydrophilic coating, metal & plastic straightner, metal trocar & drainage bag connector