

# DENSSICATH<sup>®</sup>

## MULTI-PURPOSE DRAINAGE CATHETER SET - MALECOT

Used for percutaneous placement of Malecot tip catheter in liver for biliary drainage / in the abdominal or thoracic cavity for drainage of clear or infected fluid collection / abscess.

### Set Includes

- Scalpel blade
- Initial puncture needle
- PTFE guidewire J tip
- Fascial dilators
- Malecot catheter with hydrophilic coating
- Metal & plastic straightener
- Catheter fixing disc
- Drainage bag connector



**SET**

### Size (in Fr)

8, 10, 12, 14, 16, 18, 20 & 22

### Length (in cm)

22, 25, 30 & 35

### Advantages

Kit provides entire devices for performing Seldinger type percutaneous access procedure.