

DENSSICATH[®]

MULTI-PURPOSE DRAINAGE CATHETER SET - PIGTAIL WITH LOCKING MECHANISM

Used for percutaneous placement of a pigtail tip catheter in liver for biliary drainage / in the abdominal or thoracic cavity for drainage of clear or infected fluid collection / abscess.

Set Includes

- Scalpel blade
- Initial puncture needle
- PTFE guidewire J tip
- Fascial dilators
- Pigtail catheter with locking mechanism & hydrophilic coating
- Metal & plastic straightener
- Catheter fixing disc
- Drainage bag connector



Size (in Fr)

7, 8, 9, 10, 12, 14 & 16

Length (in cm)

22, 25, 30 & 35

Advantages

- Kit provides entire devices for performing Seldinger type percutaneous access procedure
- The locking mechanism prevents accidental slippage of the catheter