

Dual Lumen Ureteral Catheter

Size (in Fr)
8, 10 & 12

Length (in cm)
50 & 70



Intended Use:

- Catheterization of the urinary tract during intraoperative procedures
- Access, advancement or exchange of ancillary devices
- Injection of contrast media

Key Features:

- Has double lumen, one lumen is for insertion over the guidewire and another lumen is for injection of contrast medium
- Optimal radiopacity for better visualization
- 1 cm marking, helps in catheter position, measurement of ureter length
- Exchange of wires without removing catheter

Advantages:

- Can inject contrast while guidewire in position
- Reduced procedural time