


DURAFLEX T STENT

REINFORCED SILICONE TUMOR STENT



Type

- Both Ends Open

Size (in Fr)

6 & 7

Length (in cm)

24 & 26

Package Includes

Stent with pusher & Guidewire (Prime[®] Wire)"

Intended Use:

Used for temporary stenting of ureter in patients with extrinsic ureteral obstructions.

Key Features:

- Manufactured by using silicone blend polyurethane for enhanced patient comfort
- Reinforced Nitinol coil layer provides superior compressive strength
- Stiff shaft wall enables stent to be placed in tight strictures
- Compatible with 0.035" guidewire