

# **DURAFLEX T STENT**

REINFORCED SILICONE TUMOR STENT

Used for temporary stenting of ureter in patients with extrinsic ureteral obstructions.

## Key Features

- Manufactured by using silicone blend polyurethane for enhanced patient comfort
- Reinforced stainless steel coil layer provides superior compressive strength
- Stiff shaft wall enables stent to be placed in tight strictures
- Compatible with 0.035" guidewire

## Also available in set, which Includes

Stent with or without string, pusher & guidewire

## Type


- Both Ends Open

## Size (in Fr)

6 & 7

## Length (in cm)

24 & 26



MRI CONDITIONAL