


DURAFLEX T STENT

REINFORCED SILICONE TUMOR STENT



MRI CONDITIONAL



Used to maintain luminal patency and flow rate during extrinsic compression of ureter

300+

Deployments
in India

Luminal Patency
Under Extrinsic
Pressure

***95%**


1+

Years Indwelling
Period


Additional Devices
Required
for Deployment

ZERO

ADDED ADVANTAGES



Metal coil ensures ample flow rate, thereby relieving obstructive uropathy in MUO (Malignant Ureteral Obstruction)



Metal coil in silicone-blended Polyurethane diminishes tissue invasion and ingrowth



Pigtail loops provide shape memory to avoid stent migration

SET INCLUDES

Duraflex T Stent, Pusher,
Prime Wire (PTFE Coated Nitinol Guidewire
with Hydrophilic Tip) and Patient Card.

BLUENEEM[®]



Blue Neem Medical Devices Pvt Ltd

Plot Nos 270 & 271, Road No 5, Harohalli Industrial Area
Phase II, Kanakapura Taluk, Ramanagara Dist - 562112
Karnataka, India



+91 80 2976 1335 / 36



marketing@blueneem.com