


## Durasil Ureteral Stent



### Types

- Both Ends Open
- One End Closed
- With or Without String

### Size (in Fr)

5, 5.5, 6, 6.5, 7 & 8

### Length (in cm)

12, 14, 16, 18, 20, 22, 24  
26, 28 & 30

### Package Includes

Stent with pusher

### Also available in set, which includes

Stent with pusher &  
guidewire

### Intended Use:

Intended to be implanted for drainage of urine from renal pelvis to urinary bladder.

### Key Features:

- Manufactured by using silicone blend polyurethane for enhanced patient comfort
- Optimal loop shape and high memory retention reduces the risk of migration
- Optimal radiopacity for better visualization