

# DURASIL URETERAL STENT

Intended to be implanted for drainage of urine from renal pelvis to urinary bladder.

## Key Features

- Manufactured by using silicone blend polyurethane for enhanced patient comfort
- Optimal loop shape and high memory retention reduces the risk of migration
- Optimal radiopacity for better visualization

## Package Includes

Stent with pusher

## Also available in set, which includes

Stent with pusher & guidewire

## Types

- Both Ends Open
- One End Closed
- With or Without String

## Size (in Fr)

5, 5.5, 6, 6.5, 7 & 8

## Length (in cm)

12, 14, 16, 18, 20, 22, 24, 26, 28 & 30

