

## Grasping Forceps (Biprong)



### Intended Use:

Used to remove the stones and retrieval of foreign objects under direct vision.

URS	
<b>Size</b> (in Fr)	3, 4 & 5
<b>Length</b> (in cm)	60

PCNL	
<b>Size</b> (in Fr)	7, 8, 9 & 12
<b>Length</b> (in cm)	40

Cystoscope	
<b>Size</b> (in Fr)	3, 4 & 5
<b>Length</b> (in cm)	40

Stent Removing Forceps	
<b>Size</b> (in Fr)	5
<b>Length</b> (in cm)	40