

## HYPOSPADIAS CALIBRATION SET

---



### Set Includes

---

Calibrators of Sizes  
5, 6, 7 & 8 Fr x 20 cm

---

PTFE Coated Guidewire  
0.025" x 40 cm

### Intended Use:

Used for gentle post-operative dilatation of the neo-urethra.

### Key Features:

- Suited For post operative calibration of Neo urethra
- Made of soft biocompatible tested Material
- Soft and atraumatic
- Easy to use set of graded sized calibrators
- Has a central lumen for use of guide wire if needed