


KIDSOFT HYPOSPADIAS STENT DIAPER (URETHRAL STENT)



Package Includes
Stent with Stylet

Variant

- Straight
- Single J

Types

- Open End
- Closed End

Size (in Fr)

6, 7, 8 & 10

Intended Use:

Used for stenting the urethra during hypospadias repair and to allow postoperative drainage of the bladder.

Key Features:

- Made of soft grade polymer for patient comfort
- Wide holes at proximal end for enhanced drainage
- Distal end provided with anchoring holes
- Pigtail at the proximal end avoids migration
- Extension tubing is provided to facilitate double diaper need