

# KIDSOFT HYPOSPADIAS STENT SINGLE J (URETHRAL STENT)



**Package Includes**  
Stent with Stylet &  
Pinch clamp

## Types

- Open End
- Closed End

**Size (in Fr)**  
6, 7, 8 & 10

### Intended Use:

Used for stenting the urethra during hypospadias repair and to allow postoperative drainage of the bladder.

### Key Features:

- Made of soft grade polymer for patient comfort
- Distal end comes with a hub that can be connected to urine bag connector
- Wide holes at proximal end for enhanced drainage
- The locking mechanism prevents accidental slippage of the stent
- Closure clamp provided to prevent urine dripping during procedure
- Self retaining single "J" at the proximal end avoids outward migration