

KIDSOFT HYPOSPADIAS STENT STRAIGHT (URETHRAL STENT)



Package Includes
Stent with Stylet &
Pinch clamp

Types

- Open End
- Closed End

Size (in Fr)

6, 7, 8 & 10

Intended Use:

Used for stenting the urethra during hypospadias repair and to allow postoperative drainage of the bladder.

Key Features:

- Made of soft grade polymer for patient comfort
- Distal end comes with trumpet shape to avoid inward migration
- Wide holes at proximal end for enhanced drainage
- Distal end provided with anchoring holes
- Closure clamp provided to prevent urine dripping during procedure