

# INITIAL PUNCTURE NEEDLE

## TWO PART

Intended to make the initial percutaneous puncture to create access and to enable the introduction of a guidewire during percutaneous interventional procedures.

### Key Features

- Siliconised needle for less trauma
- Echogenic tip for better visualization during ultrasound
- 1 cm marking on the shaft



### Type

Two Part

### Size (in Fr)

14, 16 & 18

### Length (in cm)

10, 15, 20, 23 & 25

### Advantage of Locking Hub IP Needle

During the procedure the locking hub helps in holding the inner needle with outer cannula together, which helps in easy initial puncture.