


INJECTION NEEDLE (WITH SIDE HOLES)

Used during fine needle aspiration procedure for aspiration of tissue masses, fluids etc....

Set Includes

- The special echogenic tip increases the visibility under ultra sound guidance and guarantees the needles point's localization
- Designed for percutaneous ethanol instillation to treat the small tumors
- Skinny needle with side holes and a trocar
- Side holes provide uniform instillation of fluid around the needle



Size (in G)

20

Length (in cm)

15

Advantage of Locking Hub IP Needle

The sharp point enables soft penetration