


PFlex

Ureteral Stent

Intended Use:

Used to be implanted for drainage of urine from renal pelvis to urinary bladder.



Types	Size (in Fr)	Length (in cm)
Both Ends Open & One End Closed	2.8, 3, 3.5 & 4	10, 12, 14, 16, 18 & 20

Package Includes:
Stent with pusher & guidewire

Key Features and Benefits:

- Polyurethane material enables smooth surface & low encrustation tendency
- Pigtail Loop avoids slippage & has very good retention
- Radiopaque for better visualization
- Indwelling period is 30 days