


Kidsoft[®] Hypospadias Catheter Kit Diaper

Intended Use:

Used for performing Hypospadias/Epispadias Repair & to allow post-operative drainage of urine.



Variant
Straight & Single J

Types
Open End & Closed End

Size (in Fr)
6, 7, 8 & 10

Length (in cm)
50

Set Includes:

Catheter with Pinch Cap, Flexible Stylet, Hypospadias Tourniquet, Surgical Skin Marker and Sterile Surgical Ruler

Key Features and Benefits:

- Made of soft grade polymer for patient comfort
- Wide holes at distal end for enhanced drainage
- Pigtail at the distal end avoids migration
- Extension tubing is provided to facilitate double diaper need