


Kidsoft[®] Hypospadias Catheter Kit Single J

Intended Use:

Used for performing Hypospadias/Epispadias Repair & to allow post-operative drainage of urine.



Variant
Single J

Types
Open End & Closed End

Size (in Fr)
6, 7, 8 & 10

Length (in cm)
15

Set Includes:

Catheter with Pinch Cap, Flexible Stylet, Hypospadias Tourniquet, Surgical Skin Marker and Sterile Surgical Ruler

Key Features and Benefits:

- Made of soft grade polymer for patient comfort
- Proximal end comes with trumpet shape to avoid inward migration
- Wide holes at distal end for enhanced drainage
- Proximal end provided with anchoring holes
- Pinch clamp provided to prevent urine dripping during procedure