


# Kidsoft<sup>®</sup> Hypospadias Catheter Kit Straight

## Intended Use:

Used for performing Hypospadias/Epispadias Repair & to allow post-operative drainage of urine.



**Variant**  
Straight

**Types**  
Open End & Closed End

**Size (in Fr)**  
6, 7, 8 & 10

**Length (in cm)**  
15

### Set Includes:

Catheter with Pinch Cap, Flexible Stylet, Hypospadias Tourniquet, Surgical Skin Marker and Sterile Surgical Ruler

### Key Features and Benefits:

- Made of soft grade polymer for patient comfort
- Proximal end comes with trumpet shape to avoid inward migration
- Wide holes at distal end for enhanced drainage
- Proximal end provided with anchoring holes
- Closure clamp provided to prevent urine dripping during procedure