


Kidsoft[®] Laparoscopic Pyeloplasty Stent Kit (Ureteral Stent & Set)

Intended Use:

Used for Internal Drainage of Urine from Kidney to the Bladder.
(After pyeloplasty procedure)



Size (in Fr)
2.8, 3, 3.5, 4, 5 & 6

Length (in cm)
12, 14, 15, 16, 18 & 20

Set Includes:

Introducer with Needle, PTFE Guidewire (Straight), Stent and Pusher

Key Features and Benefits:

- Ergonomically designed for stent placement during Laparoscopic Pyeloplasty
- The needle can be percutaneously introduced and placed across the anastomosis after removing the sharp stylet
- The stent can be guided into the ureter with ease, the Proximal end can be brought out as a nephrostomy if introduced through the Kidney or Pelvis
- The stent can be retained for upto 30 days