

Malecot Catheter



Package Includes

Malecot Catheter, metal or plastic Straightner & urine bag connector

Size (in Fr)

8, 10, 12, 14, 16, 18, 20 & 22

Length (in cm)

20, 25, 30 & 35

Intended Use:

Malecot Catheter is intended to be inserted percutaneously into the renal pelvis to perform external drainage of urine.

Key Features:

- The Malecot flower design is suitable for drainage of thick viscous fluids
- Integral flower design matches with the catheter size when closed and cause less trauma during insertion