

N Quadra[®] Z (Nitinol Stone Basket Zero Tip)



Size (in Fr)
2.5, 3, 4 & 5

Length (in cm)
70, 90, 110 & 120

Intended Use:

Stone Basket is intended to be used in an endoscopic procedure for removal of urinary stones or foreign objects from the Ureter or the Kidney.

Key Features:

- Zero tip design reduces chances of perforation and mucosal damage, and also avoids stone migration during procedure
- Nitinol wire construction with high strength and memory helps the basket to retain its shape during extreme torsion
- In case of stone impregnation in the basket during the procedure, ergonomic handle permits full disassembly