


## N Quadra<sup>®</sup> Flexible (Nitinol Stone Basket Flexible)

### Intended Use:

Flexible Stone Basket is intended to be used in an endoscopic procedure for removal of urinary stones or foreign objects from calyx



#### Size (in Fr)

1.8, 1.9, 2, 2.2 & 2.4

#### Length (in cm)

120

### Key Features and Benefits:

- Flexible shaft can angulate up to 270°, as per the flexible scope direction during RIRS procedure
- Kink resistant braided shaft provides durability, enabling removal of several fragments
- Tipless design helps in easy capturing of stones by allowing the basket to open at the base of calyx and encircle the stone
- Nitinol wire construction with high strength and high memory helps the basket to retain its shape during extreme torsion