

N Quadra[®] Helical Z (Nitinol Stone Basket Helical Zero Tip)

Intended Use:

Stone Basket is intended to be used in an endoscopic procedure for removal of urinary stones or foreign objects from the Ureter or the Kidney.



Size (in Fr)
2.5, 3, 4 & 5

Length (in cm)
70, 90, 110 & 120

Key Features and Benefits:

- Zero tip design reduces chances of perforation and mucosal damage, which also avoids stone migration during the procedure
- Nitinol wire construction with high strength and high memory helps the basket to retain its shape during extreme torsion
- Helical design helps in capturing small stones and fragments effectively
- In case of stone impregnation in the basket during the procedure, ergonomic handle permits full disassembly