## N QUADRA® - NITINOL STONE BASKET WITH TIP

Stone Basket is intended to be used in an endoscopic procedure for removal of urinary stones or foreign objects from the Ureter or the Kidney.

## **Key Features**

- Nitinol wire construction with high strength and high memory helps the basket to retain its shape during extreme torsion
- In case of stone impregnation in the basket during the procedure, ergonomic handle permits full disassembly
- Atraumatic tip

## Size (in Fr)

2.5, 3, 4 & 5

## Length (in cm)

70, 90, 110 & 120

