

# N QUADRA<sup>®</sup> FLEXIBLE - NITINOL STONE BASKET FLEXIBLE

Flexible Stone Basket is intended to be used in an endoscopic procedure for removal of urinary stones or foreign objects from calyx.

## Key Features

- Flexible design shaft which can angulate up to 270 degrees as per the flexible scope direction during RIRS Procedure
- Kink resistant braided shaft provides durability, making removal of several fragments possible
- Tipless design helps in easy capturing of stones by allowing the basket to open at the base of calyx and encircle the stone
- Nitinol wire construction with high strength and high memory helps the basket to retain its shape during extreme torsion
- In case of stone impregnation in the basket during the procedure, ergonomic handle permits full disassembly

## Size (in Fr)

1.8, 1.9, 2, 2.2 & 2.4

## Length (in cm)

120

