

N QUADRA[®] HELICAL Z - NITINOL STONE BASKET HELICAL ZERO TIP

Stone Basket is intended to be used in an endoscopic procedure for removal of urinary stones or foreign objects from the Ureter or the Kidney.

Key Features

- Zero tip design reduces chances of perforation and mucosal damage, which also avoids stone migration during the procedure
- Nitinol wire construction with high strength and high memory helps the basket to retain its shape during extreme torsion
- Helical design helps in capturing small stones and fragments effectively
- In case of stone impregnation in the basket during the procedure, ergonomic handle permits full disassembly

Size (in Fr)

2.5, 3, 4 & 5

Length (in cm)

70, 90, 110 & 120

