


Nephrostomy Balloon Dilator



Product Specifications

- Shaft Size: 5 Fr
- Balloon Diameter: 10 mm
- Balloon Length: 12 cm
- Catheter Length: 55 cm
- RBP: 11 atm
- Use Guidewire: 0.035 inch

Intended Use:

Used for dilation of percutaneous tract for nephrostomy procedure.

Key Features:

- Controlled inflation pressure allows for dilating force at a nephrostomy tract
- Radiopaque marker bands enhance visualization of balloon position under fluoroscopy