


Nephrostomy Catheter Set (Infant PCN Set)

Intended Use:

Used for percutaneous placement of a pigtail catheter in the renal Pelvis for nephrostomy drainage.



Size (in Fr)
6

Package Includes:

Initial Puncture Needle, PTFE Guidewire (J Tip), Dilators
Pigtail Catheter, Scalpel Blade and Urine Bag Connector

Key Features and Benefits:

- Designed for PCN in infants with small pelvis where the standard Pigtail catheters are too big to place
- The pigtail tip configuration helps in the retention of the distal tip in the renal pelvis
- Radiopaque for radiologic confirmation of the position
- Biocompatible tested material
- Measure marking on the needle and the catheter